

## The Ling Six Sound Check

The Ling Six Sound Check is a way of checking your child's hearing aids using sounds across the frequencies to ensure they are working optimally.

Using a stetoclip, listen to the hearing aid, and say each of the sounds. There should be no distortion. Alternatively, use the sounds to gauge your child's responses while they are wearing their devices. Detailed instructions on how to administer the [Lings Six-Sound Test](#) can be found here.

This check should be completed on a regular basis.

- \* The Education Authority Sensory Service is not responsible for the content of the resources linked on this page nor does it constitute endorsement

For further advice please contact the Sensory Service via phone: 028 25 661 258 or email: [sensoryservice@eani.org.uk](mailto:sensoryservice@eani.org.uk)